

Neoderm crème pour les mains Flacon-pompe de 250 ml



La crème pour les mains Neoderm est spécialement conçue pour tous les types de peau qui doivent être soignés et protégés.

Les ingrédients doux et hydratante apportent souplesse et douceur à la peau.

La crème pour les mains Neoderm n'est pas grasse au toucher et est rapidement absorbée par la peau.

- Appliquez une petite quantité de crème et massez jusqu'à ce qu'elle soit complètement absorbée par la peau
- pH neutre pour la peau

• Ingrédients : Aqua, Sorbitan Stearate, Caprylic/Capric Triglyceride, Cetyl Alcohol, Dimethicone, Phenoxyethanol, Butyrospermum Parkii Butter, Parfum, Panthenol, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Xanthan Gum, Allantoin, Glycine Soja Sterols, thylhexylglycerin, Glycine Soja Oil, Disodium PEG-5 Laurylcitrate Sulfosuccinate, Hexyl Laurate, Sodium Hydroxide, Persea Gratissima Oil, Cetearyl Alcohol, Sodium Laureth Sulfate, Biosaccharide Gum-1, Tocopherol, Ascorbyl Palmitate, Hydrogenated Palm Glycerides Citrate, Citric Acid, Sodium Benzoate